


HEARNU[®]

PRODUCT CATALOGUE



TINNI-X



FEATURES

- A Sleek and Fall-Proof Design for All-Day Comfort
- Choose From 100+ Sounds, Including Environmental, Instrumental, and Therapeutic, for Personalized Relief
- Bluetooth 5.2 Allows for Seamless Connectivity With Smartphones, Laptops, and TVs for Music, Phone Calls, and More
- Advanced Audio Notching Allows You to Easily Adjust Therapy Settings, Track Progress, and Monitor Improvements With the Tinni-X Companion Software
- Fast Charging in Just 45 Minutes and 20 Days of Standby Time
- Therapy Creation and Usage Control Dashboard for Audiologists and Professionals

APP COMPATIBLE



App Control



Audio Streaming




Rechargeable



Lightweight Design




EASY STEPS TO USE APP




1

Take a Tinnitus matching test or enter frequency from existing report




2

Create a custom therapy from different combinations of 100+ sound options



3

You can create a custom therapy or use any of default therapy



4

Play the therapy sound with your TINNI-X device while doing normal daily activities.

*Advance
Hearing Solutions
From India,
For The World*



SCAN QR TO KNOW MORE

HEAR**NU**[®]